

SILICONE GELLHORN PESSARY WITH LONG STEM

diameter 76 mm



Product Code	30056
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	2770490
CND	U089005
EAN/UPC	6972994800277

Description

Silicone Gellhorn Pessary with Long Stem - white Ø 76 mm

Used for a second to third degree uterine prolapse or procidentia. Single patient. Individually packed in a cardboard box.

Multilingual box and instructions: GB, FR, IT, ES, PT, DE, NL, SE, LV, PL, HU, SI, HR, RO, BG, GR, Arabic.

Pessary is used to support uterine prolapse in women, including cystocele and rectocele, as well as stress urinary incontinence. This medical device provides instant relief, allowing the patient to resume her normal daily activities immediately. The device must be prescribed and positioned by a qualified specialised practitioner.

Latex free. Silicone material ensures softness and flexibility.