

# SILICONE CUBE PESSARY diameter 37 mm



Product Code	30074
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	2770496
CND	U089005
EAN/UPC	6972994800932

## Description

Silicone Cube Pessary with Drain - white Ø 37 mm

Used for a third degree uterine prolapse, including procidentia, as well as a cystocele and rectocele. It should be removed nightly. Single patient. Individually packed in a cardboard box.

Multilingual box and instructions: GB, FR, IT, ES, PT, DE, NL, SE, LV, PL, HU, SI, HR, RO, BG, GR, Arabic.

Pessary is used to support uterine prolapse in women, including cystocele and rectocele, as well as stress urinary incontinence. This medical device provides instant relief, allowing the patient to resume her normal daily activities immediately. The device must be prescribed and positioned by a qualified specialised practitioner.

Latex free. Silicone material ensures softness and flexibility.