

SUTURE TRAINING KIT (Pad + instruments + sutures)



Product Code	23022
Unit of sale	1 kit
Minimum order	1
Type	No medical device
Class	
UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279230222

Description

SUTURE TRAINING KIT

Suture training kit includes: silicone pad (23020), 6 tools and 4 blades made of stainless steel, 8 sutures and a black PVC leather case.

Intended for teaching and training only, not for human or vet use.

Multilingual label and manual: GB, FR, IT, ES, PT, DE, SE, DK, GR, RO.

Kit size 19.5x22.5x3.5 cm

Standard accessories

1 Case
1 Suture pad

6 Instruments:
1 Needle holder - 12.5 cm
1 Adson forcep - 12.5 cm
1 Iris scissor - 12.5 cm
1 Scalpel handle n.3 - 12.5 cm
1 Mosquito forcep - 12.5 cm
1 Dressing forcep - 12.5 cm

4 Blades:
2 Scalpel blades n.10
2 Scalpel blades n.11

8 Sutures:
2 Nylon sutures 3/0
2 Nylon sutures 4/0
2 Silk sutures 3/0
2 Silk sutures 4/0