


ECG thermal paper 110x140mm 143s pack - white grid



Product Code	33027
Unit of sale	box of 20 pcs
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2322675
CND	Z1302
EAN/UPC	8023279330274
PARAF	988245647

Description

PHENOL FREE: completely free of phenols, it assures maximum safety to the user

Compatible with the following models of ECG:

ECG manufacturer: NIHON KOHDEN
Model: GEM 9010, 9020, 9022, 1250, 1950

Technical Specifications

Type: pack
Width: 110 mm
Length: 140 mm
Sheets: 143
Marker: yes
Grid: white