

ECG thermal paper 210x295mm x170s pack - white grid



Product Code	33028
Unit of sale	box of 1
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2322675
CND	Z1302
EAN/UPC	8023279330281

Description

PHENOL FREE: completely free of phenols, it assures maximum safety to the user

Compatible with the following models of ECG:

ECG manufacturer: NIHON KOHDEN
Model: 1500, 1550, 9110, 9130, 9132, 9320, 2250

Technical Specifications

Type: pack
Width: 210 mm
Length: 295 mm
Sheets: 170
Marker: yes
Grid: white