

DISPOSABLE STRAIGHT FORCEPS 20 cm with 3 m cable - tip 1.0 mm



Product Code	30302
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

UK-REP	Yes
CH-REP	No

RDM (NSIS)	2329323
CND	L1204
EAN/UPC	8023279303025

Description

DISPOSABLE STRAIGHT FORCEP - 20 cm - tip 1.0 mm

Disposable bipolar forceps are supplied with 3 m cable.
Valleylab plug (same connector and compatibility as 30640).
Sterilization by Ethylene Oxide (EtO).