


# SILICONE BIPOLAR CABLE 3 m with Ethicon type connector for 30310



Product Code	30315
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	No

RDM (NSIS)	2386379
CND	L1899
EAN/UPC	8023279303155

## Description

ETHICON SILICONE BIPOLAR CABLE - 3 m (10')

For bipolar scissors. Connector 2.5 mm

Compatibility: same models as 30632, 30640, 30617\*, 30643\*

\*Need adaptor 30642