


SILICONE BIPOLAR ARTERY SEALER CABLE 3 m for 30320



Product Code	30325
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	No

RDM (NSIS)	2386379
CND	L1899
EAN/UPC	8023279303254

Description

SILICONE BIPOLAR ARTERY SEALER CABLE - 3 m (10')

Connector 4 mm, Banana plug.

Compatibility: same models as 30632, 30640, 30617*, 30643* (see above).

*Need adaptor 30642