

# SET 10 ELECTRODES - 5 cm (30330-30339)



Product Code	30328
Unit of sale	kit of 10 pcs
Minimum order	1
Type	Medical device
Class	II B

UK-REP	Not required
CH-REP	Not required

RDM (NSIS)	2746097
CND	Z12010980
EAN/UPC	8033638322277
GMDN	44683

## Description

Autoclavable Electrodes Ø 2.4 mm

Reference number: SET

Length 5 cm

Autoclavable 30 times

Set of 10 electrodes

- Reference number: 1 (code 30330) - Electrode wire 0.2 - straight
- Reference number: 2 (code 30331) - Electrode wire 0.2 - angled 45°
- Reference number: 3 (code 30332) - Electrode bend - Ø 4 - straight
- Reference number: 4 (code 30333) - Electrode bend - Ø 8 - straight
- Reference number: 5 (code 30334) - Electrode hook - angled 45°
- Reference number: 6 (code 30335) - Electrode wire 1 mm - angled 45°
- Reference number: 7 (code 30336) - Electrode straight slip-knot - straight
- Reference number: 8 (code 30337) - Electrode angled slip-knot - straight
- Reference number: 9 (code 30338) - Electrode ball point - Ø 3 mm - angled 45°
- Reference number: 10 (code 30339) - Electrode ball point - Ø 3 mm - straight

Multilingual instructions: GB, FR, IT, ES PT, DE, NL, SE, LT, PL, CZ, HU, HR, RO, GR, Arabic.