

ELECTRODE BEND-STRAIGHT diam.8 mm - 5 cm



Product Code	30333
Unit of sale	box of 6 pcs
Minimum order	1
Type	Medical device
Class	II B

UK-REP	Not required
CH-REP	Not required

RDM (NSIS)	2746097
CND	Z12010980
EAN/UPC	8033638322314
GMDN	44683

Description

Autoclavable Electrodes Ø 2.4 mm

Reference number: 4
Electrode bend - Ø 8 - straight
Length 5 cm
Autoclavable 30 times

Multilingual instructions: GB, FR, IT, ES PT, DE, NL, SE, LT, PL, CZ, HU, HR, RO, GR, Arabic.