

ELECTRODE BALL POINT-ANGLED 45° - 5 cm



Product Code	30338
Unit of sale	box of 6 pcs
Minimum order	1
Type	Medical device
Class	II B

UK-REP	Not required
CH-REP	Not required

RDM (NSIS)	2746097
CND	Z12010980
EAN/UPC	8033638322376
GMDN	44683

Description

Autoclavable Electrodes Ø 2.4 mm

Reference number: 9

Electrode ball point - Ø 3 mm - angled 45°

Length 5 cm

Autoclavable 30 times

Multilingual instructions: GB, FR, IT, ES PT, DE, NL, SE, LT, PL, CZ, HU, HR, RO, GR, Arabic.