

# ELECTRODE WIRE 1mm - ANGLED 45° - 10 cm



Product Code	30345
Unit of sale	box of 6 pcs
Minimum order	1
Type	Medical device
Class	II B

UK-REP	Not required
CH-REP	Not required

RDM (NSIS)	2746097
CND	Z12010980
EAN/UPC	8033638322161
GMDN	44683

## Description

Autoclavable Electrodes Ø 2.4 mm

Reference number: 26  
Electrode wire 1 mm.- angled 45°  
Length 10 cm  
Autoclavable 30 times

Multilingual instructions: GB, FR, IT, ES PT, DE, NL, SE, LT, PL, CZ, HU, HR, RO, GR, Arabic.