

ELECTRODE BALL POINT - STRAIGHT - 10 cm



| | |
|---------------|----------------|
| Product Code | 30349 |
| Unit of sale | box of 6 pcs |
| Minimum order | 1 |
| Type | Medical device |
| Class | II B |

| | |
|--------|--------------|
| UK-REP | Not required |
| CH-REP | Not required |

| | |
|------------|---------------|
| RDM (NSIS) | 2746097 |
| CND | Z12010980 |
| EAN/UPC | 8033638322208 |
| GMDN | 44683 |

Description

Autoclavable Electrodes Ø 2.4 mm

Reference number: 30
Electrode ball point - Ø 3 mm - straight
Length 10 cm
Autoclavable 30 times

Multilingual instructions: GB, FR, IT, ES PT, DE, NL, SE, LT, PL, CZ, HU, HR, RO, GR, Arabic.