


SET 10 ELECTRODES - 10 cm (30340-30349)



Product Code	30351
Unit of sale	kit of 10 pcs
Minimum order	1
Type	Medical device
Class	II B

UK-REP	Not required
CH-REP	Not required

RDM (NSIS)	2746097
CND	Z12010980
EAN/UPC	8033638322215
GMDN	44683

Description

Autoclavable Electrodes Ø 2.4 mm

Reference number: SET
Length 10 cm
Autoclavable 30 times

- Set of 10 electrodes
- Reference number: 21 (code 30340) - Electrode wire 0.2 - straight
 - Reference number: 22 (code 30341) - Electrode wire 0.2 - angled 45°
 - Reference number: 23 (code 30342) - Electrode bend - Ø 4 - straight
 - Reference number: 24 (code 30343) - Electrode bend - Ø 8 - straight
 - Reference number: 25 (code 30344) - Electrode hook - angled 45°
 - Reference number: 26 (code 30345) - Electrode wire 1 mm - angled 45°
 - Reference number: 27 (code 30346) - Electrode straight slip-knot - straight
 - Reference number: 28 (code 30347) - Electrode angled slip-knot - straight
 - Reference number: 29 (code 30348) - Electrode ball point - Ø 3 mm - angled 45°
 - Reference number: 30 (code 30349) - Electrode ball point - Ø 3 mm - straight

Multilingual instructions: GB, FR, IT, ES PT, DE, NL, SE, LT, PL, CZ, HU, HR, RO, GR, Arabic.