

# TUNGSTEN NEEDLE ELECTRODE 5 cm - straight - disposable



Product Code	30370
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

UK-REP	Yes
CH-REP	No

RDM (NSIS)	2329251
CND	K02030101
EAN/UPC	8023279303704

## Description

Single use Sterile tungsten needle Electrodes Ø 2.4 mm

- Straight
- Length: 5 cm

Tungsten needle has extremely high melting point (> 3400°C) which helps maintain sharpness.