

TUNGSTEN NEEDLE ELECTRODE 6 cm - angled - disposable



Product Code	30373
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

UK-REP	Yes
CH-REP	No

RDM (NSIS)	2329251
CND	K02030101
EAN/UPC	8023279303735

Description

Single use Sterile tungsten needle Electrodes Ø 2.4 mm

- Angled
- Length: 6 cm

Tungsten needle has extremely high melting point (> 3400°C) which helps maintain sharpness.