

TUNGSTEN NEEDLE ELECTRODE 7 cm - straight - disposable



| | |
|---------------|----------------|
| Product Code | 30374 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | II B |

| | |
|--------|-----|
| UK-REP | Yes |
| CH-REP | No |

| | |
|------------|---------------|
| RDM (NSIS) | 2329251 |
| CND | K02030101 |
| EAN/UPC | 8023279303742 |

Description

Single use Sterile tungsten needle Electrodes Ø 2.4 mm

- Straight
- Length: 7 cm

Tungsten needle has extremely high melting point (> 3400°C) which helps maintain sharpness.