

BALLET S/S K3 ELECTROLYSIS NEEDLE



Product Code	30401
Unit of sale	box of 50 pcs
Minimum order	1
Type	No medical device
Class	

UK-REP Not required

CH-REP Not required

EAN/UPC 8023279304015

Description

BALLET STERILE DISPOSABLE ELECTROLYSIS NEEDLES.

The Ballet needle is manufactured with both a precise point and a mirror-like surface polish. The rounded tapered point of the needle will assist the electrologist in locating the follicle opening and upon insertion the needle will slide smoothly and accurately to the base of the follicle.

- K3 needle - 0.075 mm (0.003")
- Stainless steel needles: ideal for patients with strong healthy skin who respond well to electrolysis treatment.