

BALLET INSULATED K3 ELECTROLYSIS NEEDLE



Product Code	30411
Unit of sale	box of 50 pcs
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279304114

Description

BALLET STERILE DISPOSABLE ELECTROLYSIS NEEDLES.

The Ballet needle is manufactured with both a precise point and a mirror-like surface polish. The rounded tapered point of the needle will assist the electrologist in locating the follicle opening and upon insertion the needle will slide smoothly and accurately to the base of the follicle.

- K3 needle - 0.075 mm (0.003")
- Insulates needles: suggested for patients with high sensitivity to discomforts of electrolysis.