

BALL 4 mm ELECTRODE - 15 cm - autoclavable



Product Code	30446
Unit of sale	box of 5 pcs
Minimum order	1
Type	Medical device
Class	II B

RDM (NSIS)	21393
CND	L180601
EAN/UPC	8033003260937
GMDN	61876

Description

Autoclavable Electrodes Ø 2.4 mm

Reference number: 17
Ball electrode Ø 4 mm
Length 15 cm
Autoclavable 20 times

Multilingual instructions: GB, FR, IT, ES PT, DE, NL, SE, LT, PL, CZ, HU, HR, RO, GR, Arabic.