

EU NON-STICK STRAIGHT FORCEPS 20 cm - 1 mm point



| | |
|---------------|----------------|
| Product Code | 30462 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | II B |

| | |
|--------|-----|
| UK-REP | Yes |
| CH-REP | No |

| | |
|------------|---------------|
| RDM (NSIS) | 1967543 |
| CND | L1204 |
| EAN/UPC | 8964000450413 |

Description

EUROPEAN NON-STICK BIPOLAR FORCEPS

- Model: STRAIGHT
- Connector: EUROPEAN
- Length: 20 cm (7³/₄")
- Tip style: STRAIGHT
- Tip size: 1 mm

Thanks to new Ag-Tech non-stick technology, the forceps stay clean, virtually eliminating the need for cleaning during coagulation. A tip material with excellent thermal properties prevents the sticking of tissues. The non-stick effect is permanent and will not reduce with the use of forceps.