

# EU NON-STICK JEWELLER FORCEPS 11.5 cm - 0.7 mm point



Product Code	30472
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

UK-REP	Yes
CH-REP	No

RDM (NSIS)	1967540
CND	L1204
EAN/UPC	8964000450024

## Description

### EUROPEAN NON-STICK BIPOLAR FORCEPS

- Model: JEWELER
- Connector: EUROPEAN
- Length: 11.5 (4.5")
- Tip style: STRAIGHT
- Tip size: 0.7 mm

Thanks to new Ag-Tech non-stick technology, the forceps stay clean, virtually eliminating the need for cleaning during coagulation. A tip material with excellent thermal properties prevents the sticking of tissues. The non-stick effect is permanent and will not reduce with the use of forceps.