



Product Code	30474
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

UK-REP	Yes
CH-REP	No

RDM (NSIS)	2822857
CND	L1204
EAN/UPC	8964000450345

## Description

### EUROPEAN NON-STICK BIPOLAR FORCEPS

- Model: ADSON
- Connector: EUROPEAN
- Length: 12.1 (4¾")
- Tip style: STRAIGHT
- Tip size: 1 mm

Thanks to new Ag-Tech non-stick technology, the forceps stay clean, virtually eliminating the need for cleaning during coagulation. A tip material with excellent thermal properties prevents the sticking of tissues. The non-stick effect is permanent and will not reduce with the use of forceps.