

# NEEDLE ELECTRODE diam.4 mm - 5.5 cm - autoclavable



|               |                |
|---------------|----------------|
| Product Code  | 30482          |
| Unit of sale  | box of 5 pcs   |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | II B           |

|            |               |
|------------|---------------|
| RDM (NSIS) | 805750        |
| CND        | L180601       |
| EAN/UPC    | 8054242444193 |
| GMDN       | 44683         |

## Description

S/S electrodes autoclavable 100 times.  
Length 5.5 cm.