

BLADE ELECTRODE diam.4 mm - 5.5 cm - autoclavable



| | |
|---------------|----------------|
| Product Code | 30483 |
| Unit of sale | box of 5 pcs |
| Minimum order | 1 |
| Type | Medical device |
| Class | II B |

| | |
|------------|---------------|
| RDM (NSIS) | 805729 |
| CND | L180601 |
| EAN/UPC | 8054242444179 |
| GMDN | 44683 |

Description

S/S electrodes autoclavable 100 times.
Length 5.5 cm.