

# ETHICON ETHILON MONOFILAMENT SUTURES - gauge 5/0 needle 19 mm



Product Code	22305
Unit of sale	box of 36 pcs
Minimum order	1
Type	Medical device
Class	III

UK-REP	No
CH-REP	No

RDM (NSIS)	2561795
CND	H010201010101
EAN/UPC	8023279223057
GTIN	30705031099979
GMDN	38000

## Description

ETHILON - non absorbable sterile monofilament sutures

- Ethicon code: FH5171H
- Gauge: 5/0
- Metric: 1.0
- Colour: blue
- Circle: 3/8
- Needle length: 19 mm
- Needle: FS-2
- Suture length: 75 cm
- Minimum order: box of 36 pieces

Non absorbable surgical sutures are indicated for use in general soft tissue approximation and/or ligation included use in cardiovascular, ophthalmic and neurosurgical procedures.