

NEEDLE ELECTRODE - 7 cm - autoclavable



| | |
|---------------|----------------|
| Product Code | 30511 |
| Unit of sale | 1 pc |
| Minimum order | 10 |
| Type | Medical device |
| Class | II B |

| | |
|--------|-----|
| UK-REP | Yes |
| CH-REP | Yes |

| | |
|------------|---------------|
| RDM (NSIS) | 771104 |
| CND | L180601 |
| EAN/UPC | 8033003261033 |
| GMDN | 61876 |

Description

Autoclavable Electrodes Ø 2.4 mm

Reference number: 15
Needle electrode
Length 7 cm
Autoclavable 20 times

Multilingual instructions: GB, FR, IT, ES PT, DE, NL, SE, LT, PL, CZ, HU, HR, RO, GR, Arabic.