

SURGEON 1 PIN AUTOCLAVABLE HANDLE for 50, 80, 106, 108



Product Code	30515
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

UK-REP	Yes
CH-REP	No

RDM (NSIS)	2745592
CND	Z12010980
EAN/UPC	8033638321232
GMDN	44677

Description

1-PIN MONOPOLAR HANDLES FOR ELECTRODES Ø 2.4 mm

1 pin plug HANDLES for diatermo 50, 80, 106, 108

These handles are provided with a 1 Pin socket (4 mm) fitting all Gima Diatermo and most of electrosurgical units in the market. These handles do not have push-buttons.

- SURGEON 1 PIN AUTOCLAVABLE HANDLE - 100 times
Autoclavable, waterproof handle with ergonomic grip.
Provided with 3 m cable and blade electrode