



Product Code	30516
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

UK-REP	Yes
CH-REP	No

RDM (NSIS)	100185
CND	Z12010902
EAN/UPC	8023279305166
PARAF	903150858
GMDN	44776

Description

DIATERMO 106 - 50W - monopolar

This HF electrosurgery unit is particularly suitable for small operations and is recommended for dentistry, first aid and dermatology work.

It offers precision of cut and immediate coagulation and give excellent results in a wide range of applications: from uterine coagulation to microcoagulation in dermatology.

The modes of operation, CUT, COAG or BLEND, can be easily selected by a switch on the front panel.

The complete isolation between patient and the unit and strict respect of international norms guarantee total safety.

Made in Italy.

Technical Specifications

- Output power:
 - CUT 50 W/400 ohm
 - COAG 40 W/400 ohm
 - BLEND 45 W/400 ohm
- Working frequency: 600 KHz +/- 10%
- Size: 180 x 180 x h 100 mm
- Weight: 2.8 kg
- Operating Voltage: 115/230 V - 50/60 Hz switchable
- Electric input power: 90 VA
- Fuse: 2 x 1 AT
- Alarm for overheating
- Footswitch operating
- Norms: EN60601-1 Class I CF, EN 60601-1-2 93/42/EEC (Class IIb)

Standard accessories

- KIRKY single-use handle without buttons
- MB waterproof foot-switch (30577)
- Flat neutral plate with cable (30564)
- Kit of 6 short autoclavable electrodes, Ø 2.4 mm, length 5 cm
- Power cable 2 m - Schuko plug
- User manual: available in GB, IT, ES, FR