

# SINGLE USE HANDLE for 50, 80, 106, 108



Product Code	30517
Unit of sale	box of 50 pcs
Minimum order	1
Type	Medical device
Class	II B

UK-REP	Yes
RDM (NSIS)	756587
CND	K0201010501
EAN/UPC	8033003263921
GMDN	44678

## Description

1-PIN MONOPOLAR HANDLES FOR ELECTRODES Ø 2.4 mm

1 pin plug HANDLES for diatermo 50, 80, 106, 108

These handles are provided with a 1 Pin socket (4 mm) fitting all Gima Diatermo and most of electrosurgical units in the market. These handles do not have push-buttons.

• KIRKY SINGLE-USE 50-80-106-108 HANDLE - STERILE - box of 50

Single-use handle without button with blade electrode with 4 mm PIN connection.

Multilingual instructions: GB, FR, IT, ES, PT, DE, NL, PL, RO, HU, DK, SE, FI, EE, LV, RU, SK, CZ, HR, SI, GR, BG, TR, Arabic.