

AUTOCLAVABLE HANDLE for 50, 80, 106, 108



Product Code	30519
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

RDM (NSIS)	1752533
CND	Z12010980
EAN/UPC	8033003264713
GMDN	42553

Description

1-PIN MONOPOLAR HANDLES FOR ELECTRODES Ø 2.4 mm

1 pin plug HANDLES for diatermo 50, 80, 106, 108

These handles are provided with a 1 Pin socket (4 mm) fitting all Gima Diatermo and most of electrosurgical units in the market. These handles do not have push-buttons.

• AUTOCLAVABLE 50-80-106-108 HANDLE - 100 times

Autoclavable handle without buttons, provided with a 3 meter long cable and a blade electrode.