

# ORION PLUS AUTOCLAVABLE HANDLE - 100 times with screw clamp



Product Code	30539
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	1120032
CND	Z12010980
EAN/UPC	8033003262665
GMDN	42552

## Description

Autoclavable handle with 2 large "click" push-button, provided with a 3 meters long silicone cable and a blade electrode. It can be autoclaved up to 100 times (20 times electrode). For electrodes Ø 2.4 mm.

Equipped with a screw clamp locking system that not only guarantees a safe contact, but also adapts to the different electrode diameters. High insulation to prevent liquid damage.