

MB AUTOCLAVABLE HANDLE 100 times - 5 m cable



Product Code	30545
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

RDM (NSIS)	773562
CND	Z12010980
EAN/UPC	8033003261668
GMDN	42552

Description

Autoclavable handle with 2 large "click" push-button, provided with a 5 meter long silicone cable and a blade electrode. It can be autoclaved up to 100 times (20 times electrode). For electrodes Ø 2.4 mm.