

CONNECTION CABLE FOR PADS (30490-30495) - Valleylab type - 5 m



Product Code	30547
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes
RDM (NSIS)	2224052
CND	Z12010980
EAN/UPC	8033003263174
GMDN	35042

Description

CABLES TO CONNECT GROUND PADS TO ELECTROSURGERY UNITS

- CONNECTION CABLE FOR PLATES VALLEYLAB TYPE - 5 m - Universal Valleylab connection

Universal cable to connect plates (30490-30495) to Gima diatermo.
Safety lock with universal connection.

The bipolar universal connection cables are supplied with a bipolar clamp for connection to disposable ground pads with a tab (dimensions 22-25 mm).

Various connectors make the cables compatible with numerous electrosurgical devices.

The cable can be fixed to the surgical sheet during the operation thanks to a metal sliding clamp.

Latex free.

Individually packed in sealed polyethylene pouch compliance to 93/42/EEC Medical Device Class I.

COMPATIBLE ELECTROSURGERY UNITS:

- Aspenlabs: Excalibur
- Valleylab: REM type devices
- Bowa
- Conmed
- Erbe ICC International
- Erbe ACC International
- Erbe VIO International