

STAR AUTOCLAVABLE MB HANDLE 30 times - 3 m cable



Product Code	30548
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

RDM (NSIS)	21596
CND	Z12010980
EAN/UPC	8033003261606
GMDN	42552

Description

Autoclavable handle with 2 large "click" push-button, provided with a 3 meters long silicone cable and a blade electrode. It can be autoclaved up to 30 times (20 times electrode). For electrodes Ø 2.4 mm.