



Product Code	30549
Unit of sale	box of 50 pcs
Minimum order	1
Type	Medical device
Class	II B

UK-REP	Yes
--------	-----

RDM (NSIS)	513555
CND	K0201010501
EAN/UPC	8033003264188
GMDN	44681

Description

3-PIN MONOPOLAR HANDLES FOR ELECTRODES Ø 2.4 mm

- KIRKY SINGLE-USE HANDLE - STERILE - box of 50
Single-use handle with "click" push-buttons for cut and coag.
Includes blade electrode and 3.20 m cable

CAUTERY & ELECTROSURGERY 230 VALLEYLAB CONNECTOR HANDLES with 2 "click" push-buttons

GIMA electrocautery pencils with CUT-COAG switch offer maximum reliability and functionality to the operator. Provided with a 3 Pin plug Valleylab type fitting all Gima mono-bipolar Diatermo MB (80, 120, 122, 132, 160, 200, 240, 300, 380, 400) and most of bipolar electrocautery units: Aspen Lab, ES Generator, EMC Statome, Danieli, Braun, Berchtold, Ohmeda, Ritter, Conmed, Eschmann, Bovie, Martin, Kontron, Bard, Birtcher 771, Cobi, Valleylab, Led, Erbe...

Multilingual instructions: GB, FR, IT, ES, PT, DE, NL, PL, RO, HU, DK, SE, FI, EE, LV, RU, SK, CZ, HR, SI, GR, BG, TR, Arabic.