



Product Code	30553
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

UK-REP	Yes
CH-REP	No
RDM (NSIS)	511977
CND	K0201010501
EAN/UPC	8964000207017

## Description

3-PIN MONOPOLAR HANDLES FOR ELECTRODES Ø 2.4 mm

- GIMA SINGLE-USE HANDLE - STERILE  
Handle with "click" push-buttons for cut and coag.  
With 3 m cable.

CAUTERY & ELECTROSURGERY 230 VALLEYLAB CONNECTOR HANDLES with 2 "click" push-buttons  
GIMA electrosurgery pencils with CUT-COAG switch offer maximum reliability and functionality to the operator. Provided with a 3 Pin plug Valleylab type fitting all Gima mono-bipolar Diatermo MB (80, 120, 122, 132, 160, 200, 240, 300, 380, 400) and most of bipolar electrosurgical units: Aspen Lab, ES Generator, EMC Statome, Danieli, Braun, Berchtold, Ohmeda, Ritter, Conmed, Eschmann, Bovie, Martin, Kontron, Bard, Birtcher 771, Cobi, Valleylab, Led, Erbe...