

MB AUTOCLAVABLE HANDLE 100 times - 3 m cable



Product Code	30554
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

RDM (NSIS)	1729229
CND	L180301
EAN/UPC	8033695725684
PARAF	964968073
GMDN	61870

Description

Autoclavable handle with 2 large "click" push-button, provided with a 3 meters long silicone cable and a blade electrode. It can be autoclaved up to 100 times (20 times electrode). For electrodes Ø 2.4 mm.

GIMA electrocautery pencils with CUT-COAG switch have been studied and tested to offer maximum reliability and functionality to the operator. MB handles are provided with a 3 Pins socket Valley Lab type that fits all Gima mono-bipolar Diatermo MB (122, 132, 160, 200, 202, 240, 380) and majority of bipolar electrocautery units: Aspen Lab, ES Generator, EMC Statome, Danieli, Braun, Berchtold, Ohmeda, Ritter, Conmed, Eschman, Bovie, Martin, Kontron, Bard, Birtcher 771, Cobi, Valley Lab, Led, Erbe...