

# CONNECTION CABLE FOR PADS (30490-30495) - plug 6.3 mm - 5 m



Product Code	30563
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2224048
CND	Z12010980
EAN/UPC	8033003263167
GMDN	46906

## Description

### CABLES TO CONNECT GROUND PADS TO ELECTROSURGERY UNITS

- CONNECTION CABLE FOR PLATES - 5 m - Universal jack mono plug 6.3 mm

The bipolar universal connection cables are supplied with a bipolar clamp for connection to disposable ground pads with a tab (dimensions 22-25 mm).

Various connectors make the cables compatible with numerous electrosurgical devices.

The cable can be fixed to the surgical sheet during the operation thanks to a metal sliding clamp.

Latex free.

Individually packed in sealed polyethylene pouch compliance to 93/42/EEC Medical Device Class I.

### COMPATIBLE ELECTROSURGERY UNITS:

- Aesculap: TM, TM401
- Alexandra: D205, E300, 500D, 102/201, 2020/E, 202, MICRO180
- ALSA: Tutte MF100, MF360/A360, 380100, 180
- Berchtold: Elektrotom 200, 390, 400, 540
- Birtcher: 772, 773, 737, 737XL, 774
- Bovie: Bovie CSV, URO, BS, 400/400, 400SR, X10, X10L
- Burdick: U-7, SU-8
- Cameron-Miller: 26-1290
- Cobe: TS-403
- Concept: 9600, 9700, 9800, 9900
- Dolley: COBII 203T, 502, 503T, 505T
- Ems Davol: System 2000, 5000
- Faset: Serie COBI
- Gima: all models
- GPS : all models
- I Cuatro: Megeri1
- Lamidey Series Subtilec: all models
- Lamidey Series Surgilec: all models
- Led: all models
- Martin: all models
- Neomed: Serie 3000, 3000A, 3002, 3006; 3010, 3020
- Philips: TS-403
- Ritter/Mdt: Solide State 400, 400B, 400SR
- Sidevan: Seb2
- Soxil: BP 2000, X10
- Thackray: T300S
- Torry: Lem 300 CST