

BOVIE HIGH TEMPERATURE FLEXIBLE LOOP TIP - sterile



Product Code	30569
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	II B

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	1909612
CND	K0299
EAN/UPC	8023279305692

Description

INTERCHANGEABLE STERILE TIPS - 10 tips
For cauteries code 30590/30584

- High temp 2" flex loop tip - sterile
- Recommended for: GYN/ENT