

WATERPROOF FOOTSWITCH for 50,80,MB120,122,132,160,200,202 - 3 m



Product Code	30578
--------------	-------

Unit of sale	1 pc
--------------	------

Minimum order	1
---------------	---

Type	No medical device
------	-------------------

Class	
-------	--

UK-REP	No
--------	----

CH-REP	No
--------	----

EAN/UPC	8033638321607
---------	---------------

Description

MB waterproof footswitch - 3m cable

Footswitche for 50, 80, 106, MB 80D, 120D, 120F, 122, 132, 160, 160D, 200, 200D, 200F, 200 T, 202