



| | |
|---------------|----------------|
| Product Code | 30579 |
| Unit of sale | box of 10 pcs |
| Minimum order | 1 |
| Type | Medical device |
| Class | II B |

| | |
|------------|---------------|
| RDM (NSIS) | 2744405 |
| CND | Z12010901 |
| EAN/UPC | 0088748217716 |

Description

BOVIE CAUTERIES

Cautery with fixed tip

- High temp fine tip
- Recommended for: DER/Cardio

- battery operated cautery offered with changeable or fixed tip
- modified design enables clinicians to "snap off" the cautery tip for safe and effective removal of the battery
- recessed power button reduces the potential for inadvertent activation
- interchangeable sterile tips (non sterile tips for 30584 and 30590) or fixed tips for 30579, 30588
- reusable handles available (30590)
- no sterilization necessary