

BOVIE CHANGE-A-TIP HIGH TEMPERATURE CAUTERY KIT



Product Code	30584
Unit of sale	1 kit
Minimum order	1
Type	Medical device
Class	II B

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	1911546
CND	K0299
EAN/UPC	0607151045012

Description

High temp cautery kit (1 handle + 4 non sterile tips)

Change-A-Tip high-temp cautery, kit including 1 handle, 4 tips (2 as 30568 and 2 as 30571 but not sterile) and 2 AA alkaline batteries in a foam-lined case