

# CRYOMEGA CRYOSURGICAL DEVICE with 16 g cartridge



Product Code	30586
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
CH-REP	Yes
EAN/UPC	0850077006005
GMDN	44724

## Description

CRYOMEGA® WITH 16 G CARTRIDGE - non-refillable

CryOmega® dispenses liquid nitrous oxide (-89°C) the coldest portable cryogen and it is cleared for 16 common lesions including verrucae, skin tags, and seborrheic keratoses.

Ideal choice in portable cryosurgery for physician offices looking for convenience and great first-time outcomes.

It provides clinically effective freezing with pinpoint precision.

No contact device for use on different patients at 1 cm distance.

Up to 300 seconds treatment time per device (approx 30 treatments of 10 seconds).

For medical healthcare professionals.

Multilingual pack and instructions (GB, FR, IT, ES, DE, NL, PL, FI).

Made in USA.

## Technical Specifications

Suggested freezing time as reported in literature

- Verruca Vulgaris: 10-20 seconds freezing time
- Verruca Plantaris: 10-20 seconds freezing time
- Molluscum Contagiosum: 3-10 seconds freezing time
- Skin Tags: 5-10 seconds freezing time
- Lentigo: 2-5 seconds freezing time
- Actinic Keratosis: 5 seconds freezing time
- Verruca Plana: 5-15 seconds freezing time
- Seborrheic Keratosis: 5-10 seconds freezing time
- Genital Lesions: 5-12 seconds freezing time