

BOVIE HIGH TEMP CAUTERY WITH LOOP TIP - fixed tip



Product Code	30588
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	II B

RDM (NSIS)	2744405
CND	Z12010901
EAN/UPC	0088748217721

Description

BOVIE CAUTERIES

Cautery with fixed tip

- High temp loop tip
- Recommended for: DER/ENT

- battery operated cautery offered with changeable or fixed tip
- modified design enables clinicians to "snap off" the cautery tip for safe and effective removal of the battery
- recessed power button reduces the potential for inadvertent activation
- interchangeable sterile tips (non sterile tips for 30584 and 30590) or fixed tips for 30579, 30588
- reusable handles available (30590)
- no sterilization necessary