

HISTOFREEZER MIX - 2 bottles 80 ml + 24 applicators 2 mm, 36 applicators 5mm



Product Code	30589
Unit of sale	1 kit
Minimum order	1
Type	Medical device
Class	II A
UK-REP	Yes
CH-REP	Yes
RDM (NSIS)	300270
CND	M9001
EAN/UPC	0850077006173
PARAF	976340238
GMDN	11067

Description

The Histofreezer® portable cryosurgical system is a low cost, highly effective cryosurgical agent that provides safe and effective treatment for warts and other common benign skin lesions.

FOR PROFESSIONAL USE.

This therapeutic treatment makes use of local freezing for the controlled destruction or removal of living tissue. Suitable for adults and children.

Effective use on 9 benign lesion types and percentage of resolution:

- Verruca Vulgaris (Common Warts): 94%
- Verruca Plantaris (Plantar Warts): 86%
- Molluscum Contagiosum: 100%
- Lentigo (Age Spots): 100% facial / 99% non facial
- Condyloma (Genital Warts): 89%
- Seborrheic Keratosis: 95%
- Acrochordon (Skin Tags): 100%
- Verruca Plana (Flat Warts): 93%
- Actinic keratosis (Facial and Non-Facial): 91% facial / 87% non facial

Saturated 2 mm applicators will stay cold for 90 seconds.

Saturated 5 mm applicators will stay cold for 120 seconds.

This kit includes two 80 ml aerosol bottles, 24 applicators 2 mm and 36 applicators 5 mm.

Supplied with box and manual in: GB, FR, IT, DE, ES, NL, PL.

From web you can download: GR, PT, TH, Arab.