



Product Code	30592
Unit of sale	1 kit
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	293533
CND	M9001
EAN/UPC	0850077006159
PARAF	976340289
GMDN	11067

Description

The Histofreezer® portable cryosurgical system is a low cost, highly effective cryosurgical agent that provides safe and effective treatment for warts and other common benign skin lesions.

FOR PROFESSIONAL USE.

This therapeutic treatment makes use of local freezing for the controlled destruction or removal of living tissue. Suitable for adults and children.

Effective use on 9 benign lesion types and percentage of resolution:

- Verruca Vulgaris (Common Warts): 94%
- Verruca Plantaris (Plantar Warts): 86%
- Molluscum Contagiosum: 100%
- Lentigo (Age Spots): 100% facial / 99% non facial
- Condyloma(Genital Warts): 89%
- Seborrhic Keratosis: 95%
- Acrochordon (Skin Tags): 100%
- Verruca Plana (Flat Warts): 93%
- Actinic keratosis (Facial and Non-Facial): 91% facial / 87% non facial

Saturated 5 mm applicators will stay cold for 120 seconds.

This kit includes two 80 ml aerosol bottles and 52 applicators 5 mm.

Supplied with box and manual in: GB, FR, IT, DE, ES, NL, PL.

From web you can download: GR, PT, TH, Arab.