

CAUTERY HANDLE "AA"



Product Code	30598
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	14968
CND	Z12010980
EAN/UPC	8033695723376
GMDN	45631

Description

CAUTERIES - 1,200°C (2,200°F)

Cautery handle "AA" (tip not included)

- interchangeable tips with nickel-chrome wire
- replaceable batteries
- reusable handles

Wide application in gynaecology, orthopedy, ear/nose/throat, dermatology, plastic and reconstruction surgery.

Length 18 cm, weight 65 g

Made in Italy.