

# L-HOOK LAPAROSCOPIC ELECTRODE - 36 cm



|              |       |
|--------------|-------|
| Product Code | 30605 |
|--------------|-------|

|              |      |
|--------------|------|
| Unit of sale | 1 pc |
|--------------|------|

|               |   |
|---------------|---|
| Minimum order | 1 |
|---------------|---|

|      |                |
|------|----------------|
| Type | Medical device |
|------|----------------|

|       |      |
|-------|------|
| Class | II B |
|-------|------|

|        |     |
|--------|-----|
| UK-REP | Yes |
|--------|-----|

|        |    |
|--------|----|
| CH-REP | No |
|--------|----|

|            |         |
|------------|---------|
| RDM (NSIS) | 1710953 |
|------------|---------|

|     |           |
|-----|-----------|
| CND | K02030101 |
|-----|-----------|

|         |               |
|---------|---------------|
| EAN/UPC | 8964000452592 |
|---------|---------------|

## Description

Laparoscopic electrode without suction for use with Diatermo MB 240, 250, 300, 380, 400 with connecting cable 4 mm (30609).

- L-Hook electrode - 36 cm