

J-HOOK LAPAROSCOPIC ELECTRODE - 36 cm



Product Code	30606
--------------	-------

Unit of sale	1 pc
--------------	------

Minimum order	1
---------------	---

Type	Medical device
------	----------------

Class	II B
-------	------

UK-REP	Yes
--------	-----

CH-REP	No
--------	----

RDM (NSIS)	1710953
------------	---------

CND	K02030101
-----	-----------

EAN/UPC	8964000452608
---------	---------------

Description

Laparoscopic electrode without suction for use with Diatermo MB 240/250/300/380/400 with connecting cable 4 mm (30609)

- J-Hook electrode - 36 cm